



BT50 and BT30

Technical Note: T13-001

BT50 and BT30 chemical resistance chart at 20°C

■ Introduction

Liquids dispensed with the BT50 and BT30 digital burettes will be in contact, constantly, with the following materials: Borosilicate glass, (BSG), PTFE, PVDF, FEP, Alumina Ceramic (AC) and Platinum Iridium (PI). The following table is a guide to help with the queries regarding liquid compatibility.

Please note that these tables are just a guide. We recommend that if there is a question regarding liquid compatibility you should exercise caution in use and refer to other chemical tables where available. Good laboratory practice would be to rinse out the liquid handling unit at the end of each day with distilled water to prevent corrosive liquids being left in contact with the parts for too long.

■ Chemical compatibility

| CHEMICAL | BSG | PTFE | PVDF | FEP | AC | PI |
|----------------------------|-----|------------|------|-----|-----|-----------------|
| Acids | | | | | | |
| Acetic, glacial | R | | VR | | R | R |
| Acetic, 25% | R | R | R | R | R | R |
| Hydrochloric, concentrated | R | | R | | VR* | R ¹ |
| Hydrochloric, 25% | R | R | R | R | R* | R |
| Sulphuric, concentrated | R | | R | | VR* | VR ² |
| Sulphuric, 25% | R | R | R | R | R* | R |
| Nitric, concentrated | R | | VR | | VR* | VR ³ |
| Nitric, 25% | R | | R | | R | R |
| Phosphoric, 25% | R | R | R | R | R | VR ⁴ |
| Formic, 25% | R | R | R | R | R | |
| Trichloroacetic, 10% | R | R | NR | R | R | |
| Formic acid, 85% | R | R | R | R | R | |
| Arsenic acid | R | | R | | R | No data |
| Boric acid, 10% | R | R | R | R | R | No data |
| Chromic acid, 20% | R | R | R | R | R | No data |
| Hydrofluoric acid, 35% | NR | Exceptions | R | R | NR | R |
| Phosphoric acid 85% | R | R | R | R | NR | VR ⁴ |
| Nitric acid, 50% | R | R | R | R | R | R |
| Sulphuric acid, 95% | R | R | R | R | R | VR ² |
| Alkalis | | | | | | |
| Ammonium hydroxide, 25% | R | R | VR | R | R | R |
| Sodium hydroxide | R | R | NR | R | R | R |
| Potassium hydroxide | R | R | NR | R | VR | VR ⁵ |
| Sodium hydroxide | R | R | NR | R | VR | VR ⁶ |
| Alcohols | | | | | | |
| Methanol, 98% | R | R | R | | R | R |
| Ethanol, 98% | R | | R | | R | R |
| Ethanol, 70% | R | | R | | R | R |
| Isopropanol, n-Propanol | R | | R | | R | R |
| Amyl alcohol, Butanol | R | | R | | R | |
| Benzyl alcohol | R | R | R | R | R | |
| Ethylene glycol | R | R | R | R | R | |
| Propylene glycol | R | R | R | R | R | |
| Glycerol | R | R | R | R | R | |



| Hydrocarbons | | | | | | |
|---------------------------------|---|---|----|---|----|----------------|
| Hexane, Xylene | R | R | R | R | R | |
| Toluene, Benzene | R | R | R | R | R | |
| Kerosene, Gasoline | R | | R | | R | |
| Tetralin, Decalin | R | | R | | R | |
| Halogenated Hydrocarbons | | | | | | |
| Methyl chloride | R | | R | | R | |
| Chloroform | R | R | R | R | R | |
| Trichloroethylene | R | R | R | R | R | R |
| Monochlorobenzene, Freon | R | | R | | R | |
| Carbon tetrachloride | R | R | R | R | R | |
| Ketones | | | | | | |
| Acetone | R | R | NR | R | R | |
| Methyl ethyl ketone | R | R | NR | | R | |
| Isopropylacetone | R | | NR | | R | |
| Methyl isobutyl ketone | R | | VR | | R | |
| Esters | | | | | | |
| Ethyl acetate | R | R | R | | R | |
| Methyl acetate | R | | R | | R | |
| Amyl and Propyl acetate | R | | R | | R | |
| Butyl acetate | R | R | R | R | R | |
| Propylene glycol acetate | R | | R | | R | |
| 2-Ethoxyethyl acetate | R | | R | | R | |
| Methyl cellosolve acetate | R | | R | | R | |
| Benzyl benzoate | R | | R | | R | |
| Isopropyl myristate | R | | | | R | |
| Tricresyl phosphate | R | | VR | | R | |
| Oxides – Ethers | | | | | | |
| Ethyl ether | R | | R | | R | |
| 1,4 Dioxane and Tetrahydrofuran | R | R | VR | R | R | |
| Dimethylsulphoxide (DMSO) | R | R | NR | R | R | |
| Isopropyl ether | R | | R | | R | |
| Solvents with Nitrogen | | | | | | |
| Dimethyl formamide | R | R | NR | R | R | |
| Diethylacetamide | R | R | NR | | R | |
| Triethanolamine | R | | VR | | R | |
| Aniline | R | R | R | R | R | |
| Pyridine | R | R | VR | R | R | R |
| Miscellaneous | | | | | | |
| Phenol, aqueous, 10% | R | | R | | R | R |
| Formaldehyde solution, 30% | R | R | R | R | R | |
| Hydrogen peroxide, 30% | R | R | R | R | R | R ⁷ |
| Silicone oil and Mineral oil | R | | R | | R | |
| Pyridine | R | R | VR | R | R | R |
| Acetaldehyde | R | R | NR | R | R | |
| Ammonia, 25% ac. Sol. | R | R | NR | | R | R ⁸ |
| Ammonium | R | | NR | | R | |
| Calcium chloride aq. Sol | R | R | R | R | R | R |
| Chlorine | R | R | R | R | R | R ⁹ |
| Chlorobenzene | R | | R | | R | |
| Fluorinated hydrocarbons | R | | R | | VR | |
| Hexane | R | R | R | R | | |



| | | | | | | |
|---------------------------------|---|---|---|---|----------------|------------------|
| Iodine (tincture of) | R | R | R | | | R |
| Potassium chloride aq. Sol. | R | | R | | | R |
| Potassium permanganate aq. Sol. | R | | R | | | R ¹⁰ |
| Magnesium chloride aq. Sol. | R | | R | | | R |
| Methylene chloride | R | R | R | R | | R |
| Sodium carbonate | R | | R | | | VR ¹¹ |
| Sodium dichromate | R | R | R | R | | No data |
| Phenol, 100% | R | R | R | R | | R |
| Mercury | R | R | R | R | R ⁺ | NR ¹² |
| Silver nitrate | R | R | R | R | R | R ¹³ |
| Toluene | R | R | R | R | R | |
| Hydrogen peroxide, 30% | R | R | R | R | R | R ⁷ |
| Xylene | R | R | R | R | R | |
| Zinc chloride, 10% | R | R | R | R | R | R |
| Zinc sulphate, 10% | R | R | R | R | R | R |

Key:**R** = Resistant**VR** = Virtually resistant**SR** = Slightly resistant**NR** = Non-resistant**Exceptions** = Resistant with exceptions**Notes:**

* Depends on temperature

+ Up to 300°C

■ Notes on the resistance of Platinum-Iridium

The literature indicated that Aqua-Regia (3 parts hydrochloric acid: 1 part nitric acid) will cause slight attack to 10% iridium platinum. In practice, alloys containing more than 3% iridium show a great resistance to attack, unless they are in the form of a very fine powder. This usually involves the fusion of the alloy with zinc to increase the surface area.

Notes

1. **Hydrochloric acid** – in the presence of oxidising may cause slight attack on prolonged boiling.
2. **Sulphuric acid** – will dull the surface with prolonged heating at above 250 °C.
3. **Nitric acid (fuming)** – may dull the surface with prolonged heating.
4. **Phosphoric acid** – may dull the surface with prolonged heating.
5. **Potassium hydroxide** – the fused salt will cause slight attack.
6. **Sodium hydroxide** – the fused salt will cause slight attack.
7. **Hydrogen peroxide 30%** - in the presence of hydrochloric acid may cause slight attack on prolonged boiling.
8. **Ammonia** – heating in an ammonia atmosphere will darken and dull the surface, leading to a porous crystalline appearance.
9. **Chlorine** – in the presence of hydrochloric acid may cause slight attack on prolonged boiling.



10. **Potassium permanganate** – in the presence of hydrochloric acid may cause slight attack on prolonged boiling.
11. **Sodium carbonate** – the fused salt may cause slight attack.
12. **Mercury** – will readily attack at any temperature.
13. **Silver nitrate** – the fused salt may cause slight attack and discolour the surface.
14. **Organic compounds** – there is no data available on most of the organic compounds listed, it is unlikely they would have any detrimental effect but we can give no guarantee to this statement.