FS-Series, Entry-level Fermentor Bioreactor, 4 – 13.4 Liters

Easy-to-use, entry-level system, with built in controls for operation as a microbial fermentor as well as mammalian/animal cell culture bioreactor. This versatility, coupled with its ability to control up to three independent reactors from a single control station, makes it ideally suited to use in R&D labs, universities, teaching facilities, testing labs, and more.

Applications: Microbes, yeast, insect, plant or mammalian cells, in batch, fed-batch and continuous culture.

Optional Item:

- 1. Oxygen Enrichment Device
- 2. One external pump
- Fluorescent Lamp On/Off control Included 4 sets of T5 Fluorescent Lamp δ a external Timer on/off control (24hr).

Software:







	Jacket vessel (Dish bottom)	Single wall with blanket	Single wall (Plain bottom)	Single wall (Dish bottom)	Air lifter
Mammalian Cell Culture	٢	©	ø	Ø	Ø
Aerobic Microorganism Culture	©	©	9	ø	•
Micro-aerobic Microorganism	©	©	©	©	Ø
Anaerobic Microorganism	©	©	©	©	Ø
To culture Fragile Cell	•	•	Ø	Ø	©
Photosynthesis Culture	Ø	Ø	•	©	©
Plant Cell Culture	©	©	Ø	÷	©
Insect Cell Culture	٢	©	Ø	•	Ø

©: Excellent ♦: Well Ø: Not Recommended

Field	Products		
Medicine Products	Antibiotics; Hormones; Vaccine; Immune Modulators; Blood Proteins etc.,		
Agriculture Products	Veterinary Medicine; Biotec Fertilizer; Microbial Insecticides etc.,		
Food Additive Products	Amino Acids ; Vitamins ; Organic Acids etc.,		
Chemical & Energy Products	Ethanol; Glycerol; Methane etc.,		
For Environment Products	Digest Waste Microorganism		
Other Product	Bioleaching Of Metals; Genetically Engineered Spiders; Silk Proteins		
Biomass Products	Baker's Yeast; Mushrooms; Algae; Single-Cell Proteins etc.,		
Enzyme	Amylases; Celluloses; Proteases; Microbial Rennin; Lipases; Glucose Isomerase; Penicillin Acylase; Cholesterol Oxidase etc.,		
Metabolism Products	 Fermentation Products – Ethanol; Lacto Acid; Butanol; Acetone Growth Factors – Amino Acid; Vitamins; Citric Acid etc., Secondary Metabolites – Antibiotics; Alkaloids etc., Gene Products: Insulin; Human Growth Hormone; Interferons; Interleukins; Blood–Clotting Factors; Serum Albumin; Hepatitis B Vaccine; Rabies Vaccine And Diarrhea Vaccine For Pigs; Monoclonal Antibodies etc., 		

242

Dry heating system (Plain bottom type vessel)							
Vessel	Single wall vessel (Plain bottom)	<u>3L</u>	5L 7	7L	10L		
·		4 Devenilie ate alare / C		9.5	10.4		
		Borosilicate glass / 3	site ss for nead plate	ana ali ililings			
		10.4" Color louch so	creen interface				
	Communication port	Remote control thro	ough Ethernet (SCADA	<u>()</u>			
	Storage Program	Up to 59,994 progra	m for different kinds o	of condition			
Control unit	Data storage	Up to 10 data files					
	Data storage intertace	USB port					
	Cabinet material	ABS front panel and	iron painted housing)			
	Dimension	Footprint: 400 x 500 mm (W x D) ; Height: 740 mm					
	Rated voltage	110V or 220V ; 50/6	0 Hz				
	Inlet Gas Flow-meter	1 – 10 LPM 2 – 20 LPM					
Aeration	Impeller	Three adjustable Rushton-type impellers					
, lei unon	Sparger						
	Baffle	Removable 316L sto	inless steel baffles				
	Heating	Bottom heating plat	te (Base unit)				
	Cooling	Cooling coil (Interno	I heat exchanger)				
Temperature	Range	5°C above coolant	up to 90°C				
	Probe	Platinum RTD probe	(P† 100); 0 ~ 90°C				
	Control mode	Programmable 15 st	teps PID controller				
	Drive	Removable Top bru	shless motor				
Agitation	Speed Range	a. 30 – 1200 rpm is for general Fermentor; b. 5 – 300 rpm is for shear-sensitive cell line (optional)					
	Control mode	Programmable 15 steps PID controller					
ъH	Range / Probe	2-12 pH / Gel-filled electrode; Autoclavable					
pii	Control mode	PID					
	Range / Probe	0 – 200% / Polarogra	aphic DO sensor; Auto	oclavable			
DO	Control mode	PID; Cascade function to response to a. Stepping increasing or decreasing agitation speed b. Oxygen Enrichment Device (optional) c. To start substrate feeding program.					
	Transmitter	-/+ 2000 mV					
OKP	Probe	Probe Gel-filled electrode; Autoclavable (optional)					
Foam	Probe	316L stainless steel p	protector with insulate	ed Teflon tube; On/Of	f controller		
	Built-in 4 sets of pumps; One external pump (optional)						
Peristaltic		Precise Stepp	ping motor; min. spee	d is 1 rpm			
pump	Speed range	0 – 65 rpm					
	Control mode	Programmable 15 steps feeding control;					
Fxhaust			tainless steel conden		J		

Optional Item:

- 1. Low Speed Agitation Driver : 5 \sim 300 rpm
- 2. Oxygen Enrichment Device
- 3. One external pump
 - Specification :
 - Model # : MU-D01
 - Rpm and increment : 20 \sim 300 rpm ; 1 rpm
 - Motor Type : Brush-less motor
 - Number of roller : 1
- 4. Fluorescent Lamp control Included 4 sets of T5 Fluorescent Lamp and a external Timer on/off control (24 hr)
- 5. ORP probe : Gel-filled electrode; Autoclavable

- 1. Building main supply 110V- 60 Hz or 220V 50 Hz with electrical safety cutoff switch.
- 2. Autoclave for sterilization.
- 3. Cooling Circulating Water Bath for cooling purpose.
- Compressor Air to supply air, which comes with a reducer to control pressure to around 1 kg/cm². This air is to be dry, oil-free and filtered.
- 5. Computer for remove control.

Dry heating system (Dish bottom type vessel)							
Vessel	Single wall vessel (Dish bottom) with blanket						
		3L	5L	7L	10L	15L	20L
	Total volume(Liter)	4	7	9.3	13.4	16.5	23
	Material	Borosilicate glass / 316L ss for headplate and all fittings					
	Control panel	10.4" Color T	ouch screen Ir	nterface			
	Communication port	Remote con	trol through Etl	nernet (SCAD	4)		
	Storage Program	Up to 59,994	program for c	ifferent kinds	of condition		
Control unit	Data storage	Up to 10 dat	a files				
Connoronin	Data storage interface	USB port					
	Cabinet material	ABS front pa	nel and iron p	ainted housing	3		
	Dimension	Footprint: 400) x 500 mm (W	x D) ; Height:	740 mm		
	Rated voltage	110V or 220V	/; 50/60 Hz				
Aeration	Impeller	Two different type of impeller – a. 3 pcs of adjustable Rushton-type impeller b. 2 pcs of adjustable Pitched-blade impeller note: 1. Rushton-type impeller is for cell line that are not considered shear-sensitive 2. Pitched-blade impeller for shear-sensitive cell line 3. The standard system includes one type impeller. It is option for another type. It is available for customer-made impeller)					
	Sparger	Orifice ring					
-	Baffle	Removable 3	316L stainless s	teel baffles			
	Heating	Heating blan	ket				
	Cooling	Cooling coil	(Internal heat	exchanger); T	ne finger heat	exchanger is f	or 3L vessel.
Temperature	Range	5°C above coolant up to 60°C Platinum RTD probe (Pt 100);					
	Probe						
	Control mode	PID controlle	<u>r for one point</u>	or up to 15 st	eps programm	able control.	
Agitation	Speed Range	Two type of agitation motor : a. 30 – 1200 rpm is for general fermentor; b. 5 – 300 rpm is for shear-sensitive cell line (optional)					
	Control mode	PID controller for one point or up to 15 steps programmable control.					
pH	Range / Probe	2-12 pH / Gel-tilled electrode; Autoclavable					
•	Control mode		ontrol	O sonsori Aut			
DO	Control mode	0 - 200% / Polarographic DO sensor; Autoclavable PID; Cascade function to response to a. Stepping increasing or decreasing agitation speed b. Oxygen Enrichment Device (optional) c. To start substrate feeding program.					
ORP	Transmitter	-/+ 2000 mV					
-	Probe	Gel-filled ele	ctrode; Autocl	avable (optio	nal)		
Foam	Probe	316L stainless	steel protecto	or with insulate	ed Teflon tube;	On/Off contro	ller
		Built-in 4 se	ets of pumps; C	<u>)ne external p</u>	oump (optional)	
Peristaltic	Speed reverse	Precis	e stepping mo	otor; min. spee	a is 1 rpm		
pump	Control mode	Programmable 15 steps feeding control; Pump can be assigned for Acid, Base, Antifoam and Substrate					
Exhaust	316L Stainless steel condenser						

Optional Item:

1. Low Speed Agitation Driver : 5 ~ 300 rpm

- 2. Oxygen Enrichment Device
- 3. One external pump
- Specification :
- Model # : MU-D01
- Rpm and increment : 20 \sim 300 rpm ; 1 rpm
- Motor Type : Brush-less motor
- Number of roller : 1
- 4. ORP probe : Gel-filled electrode; Autoclavable

Utility Requirement:

- 1. Building main supply 110V- 60 Hz or 220V 50 Hz with electrical safety cutoff switch.
- 2. Autoclave for sterilization.
- 3. Cooling Circulating Water Bath for cooling purpose.
- Compressor Air to supply air, which comes with a reducer to control pressure to around 1 kg/cm². This air is to be dry, oil-free and filtered.
- 5. Computer for remove control.

244

	1	hermostat system	(Jacket vessel type	e)	
Vessel	Jacket vessel (Dish bottom) Total volume(Liter)	ressel ttom) 3L 5L 7L 10L ne(Liter) 4 7 9.3 13.4			10L 13.4
	Control panel	10.4" Color Touch sc	creen Interface		
Control unit	Communication port Storage Program Data storage Data storage interface Cabinet material Dimension	10.4 Color Touch screen Interface Remote control through Ethernet (SCADA) Up to 59,994 program for different kinds of condition Up to 10 data files USB port ABS front panel and iron painted housing Ecotorist (400 x 500 mm (W x D) : Height: 725 mm			
	Rated voltage	110V or 220V ; 50/6	0 Hz	1	
Aeration	Impeller	Two different type of impeller – a. 3 pcs of adjustable Rushton-type impeller b. 2 pcs of adjustable Pitched-blade impeller note: 1. Rushton-type impeller is for cell line that are not considered shear-sensitive 2. Pitched-blade impeller for shear-sensitive cell line 3. The standard system includes one type impeller. It is option for another type. It is available for customer-made impeller).			
	sparger Rafflo	Onlice ring	inlass staal bafflas		
Temperature	Control system Range Probe	Thermostat system; I (400 W heater/wate 5°C above coolant Platinum RTD probe	Built-in heat exchanger circulation pump); / up to 60°C (Pt 100)	er Automatic cooling wc	ater valve
	Control mode	Programmable 15 steps PID controller			
Agitation	Drive Speed Range	Removable Top brushless motor Two type of agitation motor: a. 30 – 1200 rpm is for general fermentor; b. 5 – 300 rpm is for shear-sensitive cell line (The standard system includes one type of agitation motor. It is option for another type.)			
	Control mode	Programmable 15 steps PID controller			
Hq	Range / Probe	2-12 pH / Gel-filled electrode; Autoclavable			
•	Control mode	PID			
DO	Control mode	 U = 200% / Polarographic DO sensor; Autoclavable PID; Cascade function to response to a. Stepping increasing or decreasing agitation speed b. Oxygen Enrichment Device (optional item) c. To start substrate feeding program. 			
Foam	Probe	316L stainless steel protector with insulated Teflon tube; On/Off controller			
	Built-in 4 sets of pumps; One external pump (optional)				
Peristaltic	Created array	Precise Stepp	oing motor; min. spee	d is 1 rpm	
pump	Control mode	Programmable 15 st Pump can be assign	teps feeding control; ned for Acid, Base, Ar	nti foam and Substrate	e
Exhaust	316L Stainless steel condenser				

Optional Item:

- 1. Low Speed Agitation Driver : 5 ~ 300 rpm
- 2. Oxygen Enrichment Device
- 3. One external pump
 - Specification :
 - Model # : MU-D01
 - Rpm and increment : 20 ~ 300 rpm ; 1 rpm
 - Motor Type : Brush-less motor
 - Number of roller : 1
- 4. ORP probe : Gel-filled electrode; Autoclavable

- 1. Building main supply 110V- 60 Hz or 220V 50 Hz with electrical safety cutoff switch.
- 2. Autoclave for sterilization.
- 3. Cooling Circulating Water Bath for cooling purpose.
- Compressor Air to supply air, which comes with a reducer to control pressure to around 1 kg/cm². This air is to be dry, oil-free and filtered.
- 5. Computer for remove control.

Thermostat system (Dish bottom type vessel)						
Vessel	Single wall vessel (Dish bottom)					
-		3L	5L	7L	10L	
-	Total volume(Liter)	4	7	9.3	13.4	
	Material	Borosilicate glass / 316L ss for headplate and all fittings				
	Control panel	10.4" Color Touch so	creen Interface			
	Communication port	Remote control thro	ugh Ethernet (SCADA	۹)		
	Storage Program	Up to 59,994 progra	m for different kinds o	of condition		
Control unit	Data storage	Up to 10 data files				
	Data storage interface	USB port				
-	Cabinet material	ABS front panel and	l iron painted housing			
	Dimension	Footprint: 400 x 500	<u>mm (W x D) ; Height:</u>	740 mm		
	Rated Voltage	110V or 220V ; 50/6		0.0		
Aeration	Impeller	Two different type of impeller – a. 3 pcs of adjustable Rushton-type impeller b. 2 pcs of adjustable Pitched-blade impeller note: 1. Rushton-type impeller is for cell line that are not considered shear-sensitive 2. Pitched-blade impeller for shear-sensitive cell line 3. The standard system includes one type impeller. It is option for another type. It is available for customer-made impeller).				
	Sparger	Orifice ring				
	Baffle	Removable 316L sto	inless steel baffles			
	Control system	Inermostat system; Built-in heat exchanger(400 W heater/water circulation pump) Automatic cooling water valve			er circulation pump);	
Temperature	Range	5°C above coolant				
-		Platinum RID probe	(PT 100);			
		Programmable 15 s	reps PID controller			
Agitation	Speed Range	Two type of agitation motor : a. 30 – 1200 rpm is for general Fermentor; b. 5 – 300 rpm is for shear-sensitive cell line (optional)				
	Control mode	Programmable 15 s	eps PID controller			
рH	Range / Probe	2-12 pH / Gel-filled	electrode; Autoclavo	able		
•	Control mode					
	Range / Probe	<u> </u>	aphic DO sensor; Auto on to response to			
Control mode Contr			tation speed			
ORP	Transmitter / Probe	-/+ 2000 mV / Gel-fill	ed electrode; Autocle	avable (optional)		
Foam	Probe	316L stainless steel p	protector with insulate	ed Teflon tube; On/Of	f controller	
	Four	built-in Easy tube loc	ding pumps; One ex	ternal pump (Optiona	11)	
Peristaltic		Precise Stepp	ping motor; min. spee	d is 1 rpm		
Peristalfic	Speed range	0 – 65 rpm				
F F	Control mode	Programmable 15 s Pump can be assign	teps feeding control; ned for Acid, Base, Ar	nti foam and Substrate	e	
Exhaust		3161 \$	tainless steel conden	iser		

Optional Item:

- 1. Low Speed Agitation Driver : 5 ~ 300 rpm
- 2. Oxygen Enrichment Device
- 3. One external pump
 - Specification :

246

- Model # : MU-D01 Rpm and increment : 20 ~ 300 rpm ; 1 rpm
- Motor Type : Brush-less motor
- Number of roller : 1
- 4. Fluorescent Lamp control Included 4 sets of T5 Fluorescent Lamp and a external Timer on/off control (24 hr)
- 5. ORP probe : Gel-filled electrode; Autoclavable

- 1. Building main supply 110V- 60 Hz or 220V 50 Hz with electrical safety cutoff switch.
- 2. Autoclave for sterilization.
- 3. Cooling Circulating Water Bath for cooling purpose.
- 4. Compressor Air to supply air, which comes with a reducer to control pressure to around 1 kg/cm². This air is to be dry, oil-free and filtered.
- 5. Computer for remove control.

Thermostat system (Jacket vassel type)							
Vessel	Jacket vessel (Dish bottom)						
	Туре	Jacket Single wall (No temp. control system)					
-	Max. working volume	6.6 Liter					
-	Total volume	6.6 Liter					
-	Material	Borosilicate glass / 316L ss for headplate and all fittings					
	Control panel	10.4" Color Touch screen Interface					
-	Communication port	Remote control through Ethernet (SCADA)					
-	Storage Program	Up to 59,994 program for different kinds of condition					
Control unit	Data storage	Up to 10 data files					
Coniroi unir	Data storage interface	USB port					
-	Cabinet material	ABS front panel and iron painted housing					
Dimension		Footprint: 400 x 500 mm (W x D) ; Height: 740 mm					
	Rated voltage	10V or 220V ; 50/60 Hz					
-	Inlet Gas Flow-meter	1 – 10 LPM					
Aeration	Tube	Adjustable draft tube					
	Sparger	Micro-sparger					
-	Control system	Thermostat system; Built-in heat exchanger (400 W heater/water circulation pump); Automatic cooling water valve					
Temperature	Range	5°C above coolant up to 60°C					
· • · · •	Probe	Platinum RTD probe (Pt 100)					
Control mode		Programmable 15 steps PID controller					
54	Range / Probe 2-12 pH / Gel-filled electrode; Autoclavable						
pii	Control mode One point control						
-	Range / Probe	0 – 200% / Polarographic DO sensor; Autoclavable					
DO	Control mode	PID; Cascade function to response to a. Stepping increasing or decreasing agitation speed b. Oxygen Enrichment Device (optional item) c. To start substrate feeding program.					
	Transmitter	-/+ 2000 mV					
OKP	Probe	Gel-filled electrode; Autoclavable (optional)					
Foam	Probe	316L stainless steel protector with insulated Teflon tube; On/Off controller					
	Four	built-in Easy tube loading pumps; One external pump (optional)					
Peristaltic		Precise Stepping motor; min. speed is 1 rpm					
pump	Speed range	0 – 65 rpm					
	Control mode	Programmable 15 steps feeding control; Pump can be assigned for Acid, Base, Anti foam and Substrate					
Exhaust	316L Stainless steel condenser						

Optional Item:

- 1. Low Speed Agitation Driver : 5 \sim 300 rpm
- 2. Oxygen Enrichment Device
- 3. One external pump
 - Specification :
 - Model # : MU-D01
 - Rpm and increment : 20 ~ 300 rpm ; 1 rpm
 - Motor Type : Brush-less motor
 - Number of roller : 1
- 4. ORP probe : Gel-filled electrode; Autoclavable

MRC.23.10.13

- 1. Building main supply 110V- 60 Hz or 220V 50 Hz with electrical safety cutoff switch.
- 2. Autoclave for sterilization.
- 3. Cooling Circulating Water Bath for cooling purpose.
- Compressor Air to supply air, which comes with a reducer to control pressure to around 1 kg/cm². This air is to be dry, oil-free and filtered.
- 5. Computer for remove control.