Bench top freeze dry system offer an economical solution for processing light loads of aqueous samples specially designed for use by researchers small footprint simple operation.

FDL-10N-50-Series, Bench Top Freeze Dryers

| FDL-10N-50: 50°C Lyophilizer (freeze-dryer) | | | | | | |
|---|---|-------|---------------------|------------------|--|--|
| | FDL-10N-50-BA FDL-10N-50-8M FDL-10N-50-TD | | | FDL-10N-50-TD-MM | | |
| Picture | | | | | | |
| Freeze temperature | -50°C | | | | | |
| Vacuum degree | | <15Pa | | | | |
| Freeze-dry area | 0.12 | m² | 0.07 m ² | | | |
| Tray load | 1.2 | 2L | 0.7L | | | |
| Water capture | 3kg/24h | | | | | |
| Sample tray | Ф200n | nm x4 | Φ180mm x3 | | | |
| Power supply | 220V 50Hz 800W | | | | | |
| Machine size (mm) | L380xW500xH340 | | | | | |
| Character | Big cold tray hole, with sample pre-freeze, stainless 304# cold tray, no pipe in cold tray, anti rust transparent plastic cover for dry chamber, security and direct-viewing stainless sample shelf, the distance between trays is adjustable international standard interface to pump. | | | | | |

FDL-10N-80-Series, Floor Freeze Dryers

| FDL-10N-80: 80°C Lyophilizer (freeze-dryer) | | | | | |
|---|---|-------------------------------|---------------------|------------------|--|
| | FDL-10N-80-BA FDL-10N-80-8M FDL-10N-80-TD | | | FDL-10N-80-TD-MM | |
| Picture | | | | | |
| Freeze temperature | -80°C | | | | |
| Vacuum degree | <15Pa | | | | |
| Freeze-dry area | 0.12 | 2 m ² | 0.07 m ² | | |
| Tray load | 1.2L | | 0.7L | | |
| Water capture | 3kg/24h | 4kg/24h | 3kg/24h | | |
| Sample tray | Φ200mn | Φ200mm x4 lays Φ180mm x3 lays | | | |
| Power supply | 220V 50Hz 800W 220V 50Hz 1300W | | 220V 50Hz 800W | | |
| Machine size (mm) | L500xW500xH750 | L492xW540xH800 | L492xW500xH750 | L500xW500xH750 | |
| Character | Big cold tray hole, with sample pre-freeze, stainless 304# cold tray, no pipe in cold tray, anti rust transparent plastic cover for dry chamber, security and direct-viewing stainless sample shelf, the distance between trays is adjustable international standard interface to pump. | | | | |

FDL-10N-50/80T, Bench Top Freeze Dryers

| FDL-10N-50 T type for multi-manifolds | | | | |
|---------------------------------------|--|--|--|--|
| | FDL-10N-50T FDL-10N-80T | | | |
| Picture | | | | |
| Brief | Basic components: 2L vacuum pump, 24pcs "T" shelf, 24pcs rubber valves. This type can be customized for your request | | | |
| Freeze temperature | -50°C -80°C | | | |
| Vacuum degree | <15Pa | | | |
| Water capture | 3kg/24h | | | |
| Freeze-dry area | 0.12m ² | | | |
| Power supply | 220V 50Hz 800W | | | |
| Machine size (mm) | L380xW500xH340 L500xW500xH750 | | | |
| Character | Big cold tray hole, with sample pre-freeze, stainless 304# cold tray, no pipe in cold tray, anti rust international standard interface to pump | | | |

FDL-18N-70-Series, Floor Freeze Dryers

| FDL-18N-70: -70°C Lyophilizer (freeze-dryer) | | | | | |
|--|---------------|---------------|---------------|---------------|------------------|
| | | | | | |
| | FDL-18N-70-NH | FDL-18N-70-BA | FDL-18N-70-TD | FDL-18N-70-8M | FDL-18N-70-TD-MM |
| | | | | | Taya alayyua |

| Heat | Without heating | Basic type | Top-down stoppering | Multi-manifolds (8 manifolds) | Top-down stoppering & Multi-manifolds | |
|--------------------|--------------------|------------|------------------------|----------------------------------|---|--|
| Heating function | No | Yes | | | | |
| Freeze temperature | < -70°C | | | | | |
| Vacuum degree | <13Pa | | | | | |
| Freeze-dry area | 0.27m ² | 0.18 | 0.11 | 0.18 | 0.11 | |
| Tray load capacity | 2.7L | 1.8L | 1.1L | 1.8L | 1.1L | |
| Water capture | 6kg/24h | | | | | |
| Sample tray | Φ200mm x6 | Φ200mm x4 | Ф200mm x3 | Ф200mm x4 | Ф200mm x3 | |
| Power supply | 220V 50Hz 1800W | | | | | |
| Machine size (mm) | L540xW640xH920 | | | | | |